

Feed
and Flush
to the
power of *e*

easier

effective

evolved



 **Kangaroo
ePump***

exciting

Now better enteral feeding begins

Introducing the

KANGAROO ePUMP Enteral Feeding Pump,

the latest innovation in enteral feeding. When patients require nutrition through continuous feeding, intermittent feeding or feeding and flushing, KANGAROO ePUMP delivers in one compact, easy-to-use device.



What does *e* mean to you?

e means easy

Easy step-by-step prompts, intuitive user interface and simplified set loading.

e means everybody

Appropriate for virtually all patients, from infant through geriatric years.

e means everywhere

Easily transportable, with battery pack for up to 15 hours of backup power.

e means enteral

Continuous or bolus enteral feeding and flushing capabilities, accurately and reliably, all from one pump.



with *e*



- Attitude independent
- “Smart” pump set technology
- “ATM-style” ease of operation
- Feeding history up to 72 hours
- DEHP-free pump sets
- Anti-free-flow pump sets
- Flushing capabilities

FEATURES

- 15 hour battery
- 15 psi
- Weight 2.34 lbs without pole clamp/
2.69 lbs with pole clamp
- Feed Rate- 1-400 in 1ml increments
- Flush Rate- 10-500 in 1 ml increments
- Flush Interval- 1-24 in 1 hr increments
- Multilingual
(English, Spanish, French and Portuguese)
- Multiple alarm options



Easy-to-use pump set

ACCESSORY PRODUCT ORDERING INFORMATION

Code Number	Description	Case Quantity
PUMPS		
382400	KANGAROO EPUMP	1
382491	POWER ADAPTER	1
382492	POLE CLAMP	1
F010484WT	BATTERY	1
PUMP SETS		
774655	PROXIMAL SPIKE SET	30
772055	500ML PUMP SET	30
773656	1000ML PUMP SET	30
776255	600ML RIGID CONTAINER PUMP SET	30
776250	1200ML RIGID CONTAINER PUMP SET	30
774669	PROXIMAL SPIKE FEED AND 1000ML FLUSH SET	30
773662	1000ML FEED SET W/ 1000ML FLUSH SET	30
714650	STERILE PROXIMAL SPIKE SET	30
712054	STERILE 500ML PUMP SET	30
713650	STERILE 1000ML PUMP SET	30
716154	STERILE 100ML BURETTE SET	30
RECERTIFICATION SET		
776150	100ML BURETTE RECERTIFICATION SET	5

KANGAROO ePUMP is a trademark of Sherwood Services AG.

© 2004 Tyco Healthcare Group LP. All rights reserved.

tyco
Healthcare

Kendall

15 Hampshire Street
Mansfield, MA 02048
1-800-962-9888
508-261-8000

www.kendallhq.com

For International locations please see
www.tycohealthcare.com

H-5637 10M 07/04